

Brochure

IsoHold™



Reduce
chair time

Reduce
aerosol

Retraction

Bite block

Visualization

Protection

IsoHold™ attaches to the saliva ejector and allows the operator to suction, retract, and protect. It enables an increased field of view of the oral cavity and prep site, reduces chair time, and decreases aerosol. It also functions as a bite block. A truly unique design, its patented design allows for a multi-functional device.



IsoHold™

IsoHold™

Its unique shape was designed by integrating input from hygienists and dental assistants to create an ergonomic device that not only allows the user to retract with ease but also reduces pressure on the user's hands while retracting.

Easy Setup

IsoHold™ installs in less than 5 seconds.

Office

Saves money by reducing chair time--over 30% reduction in procedure time.


Assistants are more involved with retracting and suctioning the procedure site. This increases access and visibility to the site for the practitioner while its unique design can prevent accidental tissue injury and keep the area free of fluids.

Practitioner

Reduce fatigue by allowing access and visibility during procedures.

Applications:

- General dentistry
- Endodontics
- Periodontics
- Bite blocks
- Ultrasonic scaling
- Root planing
- Prophylaxis
- Debridements
- Sealants
- Orthodontics





IsoHold™

Easy setup

Less than 5 second setup

No additional part or attachment to purchase



Non-Latex material

No known contraindications

One-time use disposable

Caution: Federal law restricts this device to sale by or on the order of a dentist. Use of this device is prohibited for any other purpose than recommended by the manufacturer.

IsoHold™

880743 Introductory package 12 /pk

880744 Standard package 24/pk

880745 Bulk package 96/pk

www.isohold.com

818 430 9288

Admin@isohold.com

*Saliva ejector not included.

Distributed by IsoHold™, 21712 Sherman Way, Canoga Park, CA 91303 818 430-9288 Made in USA